


Silicone Finger Mat



ContainMed, Inc.

Innovation, Inside and Out



sales@containmed.com




1841 Ludlow Ave
Indianapolis, In
46201



317-487-8800

Silicone Finger Mat



Dimensions: 504mm x 244mm x 20mm



Rounded corners and finished borders




Rounded base and unobstructed vent holes prevents moisture pooling, which promotes drainage and drying



Medical Grade Silicone rated for 10,000 hours of continuous steam sterilization



Stocked in Black and Light Blue



ContainMed, Inc.

Innovation, Inside and Out