


Modular Screw Caddies



ContainMed, Inc.

Innovation, Inside and Out



sales@containmed.com




1841 Ludlow Ave
Indianapolis, In
46201



317-487-8800

Modular Screw Caddies

SURGrip



Compatible with screws as small as Ø1.5mm



Silicone rated for 10,000 hours of continuous steam sterilization



Patented silicone gripping feature eliminates the need for a lid



Simple Installation - 2 Standard Screws



Mix and Match Plates with Multiple Base Configurations to create customized screw caddy combinations




Stackable Bases

Choose Individual Screw Plates

1.5mm • 2.0mm • 2.4mm • 2.7mm
3.0mm • 3.5mm • 4.0mm • 4.5mm
5.5mm • 6.5mm

Choose Plate Bases

60mm or 115mm Base
1 to 6 Plate Slots



ContainMed, Inc.

Innovation, Inside and Out