


Screw & Implant Caddies



ContainMed, Inc.

Innovation, Inside and Out



sales@containmed.com




1841 Ludlow Ave
Indianapolis, In
46201



317-487-8800

Screw & Implant Caddies



SURGrip



Patented silicone gripping feature eliminates the need for a lid



Silicone rated for 10,000 hours of continuous steam sterilization



Compatible with screws as small as Ø1.5mm and 4mm long

Conventional Screw Caddies



Folded Sheet Metal designs with hinged cover

Implant Caddies



Precision Machined to match the implant contour



ContainMed, Inc.

Innovation, Inside and Out