



# ContainMed, Inc.

Innovation Inside and Out

5391 W. 86<sup>th</sup> Street  
Indianapolis, IN 46268  
[www.containmed.com](http://www.containmed.com)  
317.487.8800  
[team@containmed.com](mailto:team@containmed.com)

## *Rongeur & Curette Cases*




*Cases • Trays • Implant Caddies*

Contact ContainMed • [team@containmed.com](mailto:team@containmed.com) • 317.487.8800




## Rongeur & Curette Cases

### Rongeur Cases



- Upright cases hold up to 16 instruments
- XL upright case shown:
  - 16.5L x 9.5W x 6.8”H
  - Holds up to 13” long rongeurs
- Constructed of anodized aluminum and silicone rated to 10,000 hours of continuous pre-vacuum steam sterilization
- Custom sizes, color anodizing, and graphics available
- Custom, lay-flat designs also available



### Curette Cases



- 10” Curette case shown:
  - 12L x 7.6W x 3”H
  - 10.5” internal clearance
  - 0.87” diameter handles
  - Holds up to 12 instruments
- Constructed of anodized aluminum and silicone rated to 10,000 hours of continuous pre-vacuum steam sterilization
- Custom sizes, color anodizing, and graphics available

